

# **CHAPTER 11**

## **Virology**

**Ray Brunson**  
**Olympia Fish Health Center**  
**Olympia, Washington**

### **Acknowledgements**

This section is based on the Alaska Department of Fish and Game, Fish Pathology Lab Manual (Meyers, 1997) that was generously offered by Dr. Ted Meyers. We have edited, updated and adapted their manual to reflect the needs of the National Wild Fish Health Survey. The USFWS recognizes and thanks the effort and Dr. Meyers, Jill Follett, Joan Thomas, Marie Fried, Sally Short, and Tamara Burton, for their research and compilation of much of the information that is in this section.

## I. Introduction

Detection of aquatic animal viruses historically has been by growth and isolation on living cell cultures appropriately researched and chosen for the propagation of target viruses and species of host. Viral detection can also include immunological and nucleotide testing procedures. The determination of a testing procedure is a complex decision involving factors of cost, timeliness, sensitivity, specificity, efficiency, and available host tissues and technology.

For the purposes of the Wild Fish Health Survey, the USFWS has chosen the use of cell culture for initial screening and corroboration of test results using appropriate nucleotide primers for the polymerase chain reaction (PCR) technique. When PCR primers are not appropriate, immunoblot, serum neutralization, or fluorescent antibody techniques are used, see Chapter 12 - Corroborative Testing of Viral Isolates. The following sections describe the procedures and methods for virology testing by tissue culture.

**Definitions:** Several terms are used routinely in virology and throughout this section. A full Glossary of terms can be found in Appendix A.

**Media Formulations:** See Appendix B: Media Used in Tissue Culture and Virology - page 11.17.

## II. Selection of Appropriate Cell Lines

All viral testing will utilize American Type Culture Collection (ATCC) cell lines tested annually for viral sensitivity and mycoplasma infections: see section VI. Quality Control in Tissue Culture, in Chapter 10 -Tissue Culture of Fish Cell Lines.

All viral testing of salmonids and herring will be on both chinook salmon embryo (CHSE-214) and *epithelioma papillosum cyprini* (EPC) cell lines at a minimum.

All viral testing on ictalurids will be on fat head minnow (FHM) or brown bullhead (BB) and chinook salmon embryo (CHSE-214) cell lines at a minimum.

All viral testing on acipenserids will be on white sturgeon skin 1 (WSSK-1) and chinook salmon embryo (CHSE-214) cell lines at a minimum.

All viral testing for other species will be on cells selected for detection of viruses important for the specific aquatic ecosystem where samples are derived. The testing laboratory will justify the choice of cell line for alternate species, following the guidelines in Table 1.

**Table 1 - Recommended cell lines to detect target viruses in various fish species**

Family Group	1st Cell Line	Incubation Temp (C)	2nd Cell Line	Incubation Temp (C)	Virus Groups <sup>a</sup>
Acipenseridae (Sturgeon)	CHSE-214	15	WSSK-1	20	IPNV, WSI, WSH
Polyodontidae (Paddlefish)	CHSE-214	15	WSSK-1	20	IPNV
Percichthyidae (Temperate Basses)	CHSE-214	15	FHM	25	IPNV, LBBV
Cyprinidae (Carp/minnows)	CHSE-214	15	FHM	25	IPNV
Catostomidae (Suckers)	CHSE-214	15	FHM	25	IPNV
Centrarchidae (Sunfishes)	CHSE-214	15	FHM/BF-2	25	IPNV, LMBV
Percidae (Perch)	CHSE-214	15	FHM	25	IPNV, LMBV
Salmonidae (Trout / Salmon)	CHSE-214	15	EPC	15	IPNV, IHN, VHS, OMV
Clupeidae (Herring)	CHSE-214	15	EPC	15	IPNV, IHN, VHS, OMV
Ictaluridae (Catfish)	CHSE-214	15	BB/FHM	25	IPNV
Other	CHSE-21	15	FHM	25	IPNV

<sup>a</sup> Viruses: IPNV- Infectious Pancreatic Necrosis Virus and other related birnaviruses; WSI - White Sturgeon Iridovirus; WSH - White Sturgeon Herpesvirus; IHN- Infectious Hematopoietic Necrosis Virus; VHSV - Viral Hemorrhagic Septicemia Virus; OMV - *Oncorhynchus masou* Virus; LMBV – Largemouth Bass Virus and other related iridoviruses.

### III. Target Tissues

For purposes of the Survey, most tissue collection will focus on organ samples from sub-adult animals. If sexually mature adults are available, the additional testing of coelomic fluid (ovarian fluid) from gravid females or seminal fluid from gravid males may also be done. Individual samples are encouraged, but samples may be pooled if necessary.

No more than 5 individuals may be pooled into a single sample and samples from individuals must be approximately the same size. The following outline summarizes suggested sampling but the individual situation will dictate the best sample for viral testing.

**A. JUVENILE FISH  $\leq$  7.0 cm IN LENGTH**

1. Fish 2.5 cm or less in length - use whole, cut off and discard any visible yolk sac
2. Fish 2.5 to 4.0 cm - cut off and discard heads and tails
3. Fish 4.0 to 7.0 cm - use viscera

**B. JUVENILE FISH  $>$  7.0 cm IN LENGTH**

1. Kidney and spleen

**C. ADULT FISH**

1. Females - kidney and spleen, ovarian fluid from spawning or post-spawning fish.
2. Males - kidney and spleen, seminal fluid can be useful if organs are unavailable.

## **IV. Tissue Collection Procedures**

**A. WHOLE ALEVIN OR JUVENILE FISHES  $<$ 7 cm IN LENGTH**

1. Place 1 to 5 fish samples as described above into each Whirl-Pak7 bag or sterile snap cap tube per pooled sample. Keep all samples cool during the collection procedure.

**B. TISSUE SAMPLES FROM FISH  $>$  7 cm IN LENGTH**

1. Aseptically, remove a piece of the kidney and spleen from each fish with forceps, scissors, scalpel and/or tongue depressor. For pooled samples, combine tissue samples from each fish into a single Whirl-Pak7 bag or snap cap tube. The total sample for each bag or tube should be at least 0.5 g of tissue.
2. Seal Whirl-Pak7 bag or tube. Keep samples cool while collecting remainder of samples.
3. Between each fish or pooled sample, clean instruments of any tissue with gauze sponges dipped in 70% ETOH, alcohol rinse, or wipe with iodophor followed by sterile water rinse.

**C. PROCEDURES FOR COLLECTING COELOMIC (OVARIAN) FLUID SAMPLES**

1. Disinfect the abdomen of the fish with iodophor and wipe with a clean paper towel to remove any disinfectant or mucous which could drip into the sample.
2. Partially strip ovarian fluid from one female fish into a clean paper cup. If possible, avoid extrusion of blood, fecal material, and nematodes.
3. Crimp edge of paper cup to "strain out" any eggs present and pour 2-5 ml ovarian fluid from each fish into one tube per fish if samples are individual. Pour 1-2 ml per per fish if samples are pooled. Do not fill tube completely full. Do not palm or warm fluid, which could inactivate low levels of virus if present.
4. Tightly cap and place in tube rack. Keep all filled tubes in cooler chest on blue or wet

ice while collecting remainder of samples.

**An Alternative Method** - Drawing ovarian fluid using an automatic pipettor if *in situ* contamination is a problem.

- a. Install sterile tip on pipettor. Insert tip into the uro-genital opening of the fish while applying light pressure to the body.
- b. Draw up a sample of ovarian fluid and place in one sterile tube per fish or pooled sample. Discard tips between samples.

#### **D. PROCEDURE FOR COLLECTING SEMINAL FLUID:**

1. Express seminal fluid from gravid male into clean paper cup.
2. Pour 2-5 ml into tube for individual samples. Pour 1-2 ml into tube for pooled samples.
3. Tightly cap and place in tube rack.
4. Keep all filled tubes in cooler chest on blue or wet ice while collecting remainder of samples.

### **V. Transport of Tissue Samples**

#### **A. PACKING AND SHIPPING SAMPLES**

1. Check seals on bags or tubes to ensure closure.
2. Place Whirl-Pak® bags into large plastic bag. Label bag with number of samples, location sample taken, sample type, date, life stage, species of fish and enclose a completed sample Submission Form (Chapter 2 – Appendix B).
3. Place tubes in proper rack, enclose in large plastic bag and label as above
4. Transport all samples on blue ice in an ice chest. Keep samples cool and avoid freezing or sunlight.
5. Check with receiving laboratory on special shipping instructions if samples are not delivered in person.

### **VI. Processing Tissue, Coelomic and Seminal Fluid Samples**

All samples for viral testing should be processed within 48 hours and inoculated onto cell lines within 72 hours of collection.

#### **A. PROCESSING TISSUE SAMPLES**

1. Tare empty tube or bag. Weigh contents of each sample container. Keep viral samples cold during processing by placing on gel ice or cold plates.

2. Add tissue culture grade PBS or MEM-O to equal a 1:10 dilution (w/v). If toxicity is likely or suspected, additional tissues dilutions from 1:20 to 1:100 can be made. Final dilution prior to inoculation of tissue samples onto cell cultures must not exceed 1:100 (v/v).
3. Homogenize samples using a Stomacher® (Virtis® or Contorque® grinders require considerable disinfection of containers between samples). Pour or pipet ~ 3-4 mls into a 12x75mm snap-cap tube.
4. Centrifuge tissue samples at 2000x g for 20 minutes.
5. Without disturbing the pellet, aseptically pipette 1 ml supernatant from each sample into its respective tube containing 1 ml of thawed antibiotic cocktail (see Appendix B - Media Used for Tissue Culture and Virology) and vortex. If there is remaining yolk from sac fry, collect the supernatant from below the floating yolk material. Label tubes.
6. Incubate at 15°C for 2-6 hr or at 4°C overnight.

#### **B. PROCESSING COELOMIC (OVARIAN) OR SEMINAL FLUID**

1. Centrifuge ovarian or seminal fluid samples at 2000x g for 20 min if using polypropylene tubes. It is optimal to process in a refrigerated centrifuge at 4°C to prevent warming of the sample. Tubes must be balanced.
2. Without disturbing the pellet, aseptically pipette 1 ml supernatant from each sample into its respective tube containing 1 ml of thawed antibiotic cocktail (see Appendix B - Media for Tissue Culture and Virology) and vortex. Label tubes.
3. Incubate at 15°C for 2-6 hr or at 4°C overnight.

## **VII. Preparing Viral Test Plates**

The quantal assay (also referred to as endpoint dilution) is used to examine fish when only the presence or absence of a virus needs to be verified which is the purpose of the Survey. Flat-bottomed 24 well plates are usually used for this assay but other cell culture plates may be used if applicable. For determination of the Tissue Culture Infective Dose - 50% endpoint (TCID<sub>50</sub>) of a virus sample or isolate, replicate samples are necessary and 96 well plates become more useful. The TCID<sub>50</sub> assay is not routinely used in the Survey because the numbers of replicate dilutions required are often not practical. Thus, no methods will be included in this manual. Reed and Muench (1938) and Rovozzo and Burke (1973) describe the procedures for the TCID<sub>50</sub> assay. The plaque assay is another quantification method that determines plaque forming units (PFU) or infectious particles (I.P.) of a sample. Flat-bottomed 24 well plates may be used for this, but several dilutions are necessary to accurately assess the titer of PFUs. The Survey has determined not to use this test for screening. Burke and Mulcahy (1980) describe the procedure for the plaque assay.

#### **A. SEEDING FLAT-BOTTOMED 24 (16mm) WELL PLATES**

1. Determine number of plates needed for the assay. Remember that samples will require at least two cell lines (i.e. EPC and CHSE-214 for salmonids). Samples should be done in duplicate on each cell line. Additionally, both monolayer and sham controls are needed, at least one set per accession number or more if several plates are used.
2. Remove confluent cell monolayer from a tissue flask using the methods described in the section on Maintenance of Stock Cell Lines – Passage of Confluent Cell Monolayers. One confluent cell monolayer in a 75-cm<sup>2</sup> flask can make six to eight 24-well plates. Normal cell density is  $5 \times 10^5$  to  $1 \times 10^6$  cells per well, but individual labs or different cell cultures may vary.
3. Pipette 0.5 – 1.0ml cell suspension into each well of the plate(s). Mix the cell suspension frequently to keep the cells homogeneously suspended.
4. Add a few extra drops of MEM-10 to all corner wells to compensate for evaporation.
5. Any liquid spilled between wells may be aspirated off or dried by use of sterile gauze. Cover each plate with the accompanying lid.
6. Label each plate with the date, cell line initials, passage number, and operator initials. Seal lid to base with tape and place plate(s) in a plastic bag or into an airtight plastic container.
7. Incubate at 21-22°C until 100% confluent without changing the medium. The plates should be confluent overnight. If necessary, plates can be made the same day as inoculation but they will have to be seeded with more cells. However, same day inoculation does not necessarily result in any earlier detection of virus.

### **VIII. Inoculating Plates with Samples**

#### **A. MATERIALS**

1. Appropriate number of 24 well plates with 100% confluent cell monolayers (do not use plates in which cells have been confluent more than 3 days). EPC cell monolayers may require thicker confluency to avoid retraction when placed at 15°C.
2. Ovarian fluid and/or tissue samples, either in antibiotic cocktail or filtered.
3. MEM-5 dilution blanks or MEM-5 with antibiotic blanks.
4. One hundred microliter pipettor, Pipetman P100, or equivalent.
5. Sterile pipette tips.
6. Plastic bag, tape and/or airtight plastic container.

## **B. INOCULATING CELLS WITH SAMPLES**

1. Each sample is inoculated in duplicate. Label each plate with the inoculation date and case number. The assay is usually done with the plate aligned with 6 columns across the top and 4 rows down. Label wells with sample numbers and identify controls.
2. Empty medium from wells by inverting the plates over a bleach bucket or use a pipet to aspirate off, leaving a small amount of medium on the cells to prevent drying in the center.
3. Pipet 0.1 ml of each ovarian fluid or tissue sample into the replicate wells allocated per sample.
4. Inoculate control wells with 0.1 ml of MEM-0 with antibiotic cocktail diluted to the same concentration as the test samples.
5. Adsorb samples for 30 – 60 minutes at 15°C. An orbital shaker plate or rocker at 30-40 rpm is recommended.
6. Add approximately 1.0-ml of MEM-10 or MEM-5 to each well. Completely seal plate lid to base with tape or seal in a plastic bag or airtight plastic container.
7. Depending upon the virus to be detected, incubate the plate(s) at the appropriate temperature.
8. Freeze all ovarian fluids and all tissue samples at –80°C until completion of the assay.

## **C. MINIMUM LEVELS OF DETECTION (assuming replicate wells)**

1. For tissues, it is 50 infectious particles (I.P.)/g pooled sample or 250 I.P./ml/fish.
2. For ovarian fluid, not pooled, it is 10 I.P./ml
3. For ovarian fluids, pooled, it is 10 I.P./ml pooled sample or 50 I.P./ml/fish for a 5 fish pool.

## **IX. Viral Plate Observation**

Following inoculation of plates, all wells will be monitored on the following day and every other day for the next 2 working weeks for signs of CPE or toxicity. Plates will be monitored twice during the following week. Total observation period will be 3 weeks (21 days).

If no CPE, toxicity or abnormalities are observed in the 21 days, the samples are discarded and recorded as negative.



**A. Re-Inoculation** - If toxicity, abnormal pH or CPE is observed, one of the replicate wells of that sample will be aseptically aspirated, diluted 1:10 with MEM-0, filtered through a 0.45µm filter and re-inoculated onto another 24-well test plate in duplicate and monitored for an additional 14 days for CPE. All observations will be documented and recorded by the observer and kept on file with the laboratory records. If no CPE is observed in 14 days after re-inoculation, the sample is discarded and recorded as negative.

## **B. Cytopathic Effects (CPE) of Virus Infection in Tissue Culture Cells**

### **IHNV-induced CPE**

1. Rounded and granular cells in grape-like clusters.
2. Margination of nuclear chromatin (optical density of nucleoli increases and nuclear membranes appear thickened).
3. Plaques in the confluent cell monolayers are ragged in outline and contain coarsely granular debris and rounded cells. Rounded, infected cells also accumulate at plaque margins and can be present within the plaque.

### **IPNV-induced CPE**

1. Spindle-shaped or "balloon-on-a-stick"-shaped cells.
2. Pyknosis of nuclei (nuclei shrink in size and chromatin condenses).
3. Plaques are stellate in a confluent cell monolayer and contain not only live cells but also normal looking cells (these are persistently infected and will reform a normal monolayer that is virus positive).
4. Little cellular debris.

### **Herpesvirus-induced CPE**

1. Pyknosis of nuclei and cellular fusion (syncytia).
2. Syncytia produce multinucleated giant cells.
3. Plaques tend to elongate and follow whorl lines of growth if on RTG-2 cells. They have relatively clear interiors, but living cells extend into the open area.
4. Little cellular debris.

### **VHSV-induced CPE**

1. The North American VHSV isolates plaque very similarly to IHNV in EPC cells forming rounded and granular cells in grape-like clusters.
2. The European VHSV isolates differ from IHNV on RTG-2 cells by having more regular plaque margins with uniformly distributed granular debris within the plaques. Also, affected cells do not show margination of chromatin.
3. Number of days following infection with virus that CPE is usually observed in freshly monolayered fish cell cultures:

### **LMBV-induced CPE**

1. CPE within 48 hours after inoculation
2. Initial CPE - few pyknotic cells, which develop to form circular, cell free areas, with rounded cells at the margins.

3. Advanced CPE - Pyknosis, rounding and detached cell sheet. Entire cell sheet affected.

<u>Virus</u>	<u>Days</u>	
IHNV, VHSV	2-10	If blind passages are needed, do at 14 days and examine for an additional 14 days.
IPNV	2-10	" "
LMBV	2-10	" "
Herpesvirus	14-30	Frequently does not appear on the initial culture and must be blind passaged at 28 days. Examine for an additional 28 days.

**Note:** Toxicity can sometimes mimic viral CPE. Observing the gradual development of plaques over several days is the best way to distinguish viral CPE from toxicity

### C. Intensity of CPE

Monolayers are examined with an inverted light microscope at low power (125X) to determine intensity of CPE. This general scoring scheme is used to record CPE intensity:

- +1 = Only one field observed contains CPE
- +2 = Two or more fields observed contain CPE
- +3 = All fields observed contain CPE
- +4 = CPE throughout entire monolayer or monolayer no longer attached to flask/plate. (Note: Separation or retraction of cell monolayer from flask/plate edge can be due to toxicity rather than viral CPE.)

### D. Corroborative Testing

Appropriate PCR or immunological methods for corroboration and identification of virus will test CPE that is observed after re-inoculation. See Chapter 12 - Corroborative Testing of Viral Isolates, for specific protocols.

## X. Storing, Freezing and Thawing Viral Isolates

### A. PREPARATION OF VIRUS ISOLATES FOR FREEZING

1. Virus samples - suspected virus isolates from all fish species are frozen after completion of viral assays. At least 2 viral isolates (if 2 or more samples produce CPE) are frozen per date, location and species, preferably from wells having 4+ CPE.
2. Aseptically pipette 1.5-2 ml of tissue culture fluid and cell debris from the wells representing each isolate into four freezer vials. Seal tightly and label.

## **B. FREEZING VIRUS ISOLATES**

1. Freeze vials at -80°C. Virus should not be frozen in the liquid nitrogen dewer that contains the stock cell lines unless a herpesvirus is strongly suspected (i.e., the virus in whole cells could be more easily lost at -80°C).
2. Label each freezer vial per isolate with the case number, isolate number, number of passages through which cell line, fish stock and species, original sample type (ovarian fluid or tissue sample) and date frozen.
3. Log each isolate in the freezer notebook.

## **C. THAWING VIRUS SAMPLES**

1. One vial of the virus should be thawed and tested for viability before freezing all samples if the identity and stability of the isolate is unknown.\*
2. Always thaw virus isolates rapidly in lukewarm water, removing the vial just before the last of the ice in the vial has melted.
3. Decant MEM-10 from the required number of 1 to 2-day-old monolayers in 25-cm<sup>2</sup> flasks. Pipet 0.1-ml virus sample onto each cell monolayer.
4. Allow virus to adsorb for 30 minutes at 15°C.
5. Add 5 ml MEM-10 to each flask and incubate at the appropriate temperature until all cells lift off each flask (4+ CPE).

\* This may not be feasible for unknown virus isolates requiring long incubation times to produce CPE.

## XI. References

Adams, R. L. P. 1980. Cell culture for biochemists. In T. S. Work and R. H. Burdon (editors.) Laboratory techniques in biochemistry and molecular biology, Volume 8. Elsevier Northland Biomedical Press, Amsterdam. 292 p.

American Type Culture Collection. 1985. Quality control methods for cell lines, first edition.

Amos, K. H., editor. 1985. Procedures for the detection and identification of certain fish pathogens, 3rd edition. Fish Health Section, American Fisheries Society, Corvallis, Oregon. 114 p.

Batts, W.N., and J.R. Winton. 1989. Concentration of infectious hematopoietic necrosis virus from water samples by tangential flow filtration and polyethylene glycol precipitation. Canadian Journal of Fisheries and Aquatic Sciences. 46:964-968.

Batts, W.N., C.K. Arakawa, J. Bernard, and J.R. Winton. 1993. Isolates of viral hemorrhagic septicemia virus from North America and Europe can be detected and distinguished by DNA probes. Diseases of Aquatic Organisms. 17:67-71.

Burke, J. A. 1982. Aspects of the pathogenesis of infectious hematopoietic necrosis virus in two-year-old sockeye salmon *Oncorhynchus nerka*. Doctoral thesis. University of Washington, Seattle, 60 p.

Burke, J. A., and D. Mulcahy. 1980. Plaquing procedure for infectious hematopoietic necrosis virus. Applied Environmental Microbiology. 39: 872-876.

Drolet, B. S., J.S. Rohovec, and J.C. Leong. 1993. Serological identification of infectious hematopoietic necrosis virus in fixed tissue culture cells by alkaline phosphatase immunocytochemistry. Journal of Aquatic Animal Health 5:265-269.

Gravell, M., and R. G. Malsberger. 1965. A permanent cell line from the fathead minnow (*Pimephales promelas*). Annals of the New York Academy of Sciences 126: 555-565.

Meyers, T. R., J. B. Thomas, J. E. Follett, and R. R. Saft. 1990. Infectious hematopoietic necrosis virus: trends in prevalence and the risk management approach in Alaskan sockeye salmon culture. Journal of Aquatic Animal Health 2:85-98.

Meyers, T.R. (ed.). 1997. Fish pathology section laboratory manual. Alaska Department of Fish and Game, Commercial Fisheries Management and Development, Special Publication No. 12. Juneau, 187 p.

Mulcahy, D. M., G. L. Tebbit, W. J. Groberg, Jr., J. S. McMichael, J. R. Winton, R. P. Hedrick, M. Philippon-Fried, K. S. Pilcher, and J. L. Fryer. 1980. The occurrence and distribution of salmonid viruses in Oregon. Oregon State University Sea Grant College Program No. ORESU-T-80-004. 71 p.

Pacific Northwest Fish Health Protection Committee. 1987. Model Comprehensive Fish Health Protection Program.

Paul, J., 1975. Cell and tissue culture, 5th edition. Churchill Livingstone, London. 484 p.

Philippon-Fried, M. 1980. Partial characterization of six established salmonid cell lines. Master's thesis, Oregon State University, Corvallis, 58 p.

Pilcher, K. S., and J. L. Fryer. 1980. The viral diseases of fish: a review through 1978. Oregon State University Sea Grant College Program No. ORESU-R-80-019. 364 p.

Plumb, J.A, J.M Grizzle, H.E. Young, A.D. Noyes, and S. Lamprecht. 1996. An iridovirus isolated from wild largemouth bass. Journal of Aquatic Animal Health. 8: 265.

Reed, L. J., and H. Muench. 1938. A simple method of estimating fifty percent endpoints. American Journal of Hygiene 27:493-497.

Rovozzo, G. C. and C. N. Burke. 1973. A Manual of Basic Virological Techniques, Prentice-Hall, New Jersey. 287 p.

Thoesen, J. C., (ed.). 1994. Suggested procedures for the detection and identification of certain finfish and shellfish pathogens (4th Edition). Fish Health Section, American Fisheries Society, Bethesda, MD.

Watanabe, R. A., J. L. Fryer, and J. S. Rohovec. 1988. Molecular filtration for recovery of waterborne viruses of fish. Applied and Environmental Microbiology. 54:1606-1609.

Wingfield, W. H., J. L. Fryer and K. S. Pilcher. 1969. Properties of the sockeye salmon virus (Oregon strain). Proceedings for the Society of Experimental Biology and Medicine 130: 1055-1059.

Wolf, K., 1988. Fish viruses and fish viral diseases. Cornell University Press, New York. 476 p.

Wolf, K., and M. C. Quimby. 1962. Established eurythermic line of fish cells *in vitro*. Science 135:1065-1066.

Wolf, K., and T. Sano. 1975. Herpesvirus disease of salmonids. U. S. Department of the Interior, Fish and Wildlife Service, Fish Disease Leaflet 44. Washington, DC. 8 p.

## **Appendix A - Glossary of Terms used in Tissue Culture & Virology**

**BCS** - bovine calf serum taken from newborn calves.

**Blind passage** - transfer of supernatant and inoculated tissue culture cells which are not demonstrating CPE to another plate containing fresh cells in order to dilute out possible inhibitors of viral expression and/or allow possible early viral replication due to low concentrations of virus particles to progress to detectable CPE.

**Closed System** - a system of incubating cells that is sealed against the transfer of air, i.e., a flask.

**Confluent Monolayer (100%)** - a single layer of tissue culture cells in which the cells have filled in all the spaces between them.

### **Controls**

- A. Monolayer control: tissue culture cells are grown in presence of growth medium MEM-10. If CPE appears in monolayer control wells, test is invalidated and must be repeated.
- B. Sham control: diluent (MEM-0) used for suspension of samples or dilution blanks is added to cells. After adsorption, an overlay medium or MEM-10 is added. If CPE appears in sham control wells, test is invalidated and must be repeated.

**Cytopathic Effects (CPE)** - changes in the morphology and metabolism of tissue culture cells due to suspected viral infection.

**Defective Interfering Particles (D.I. Particles)** - defective or incomplete virus particles which cannot replicate but may prevent expression of the infectious virus by attaching to the tissue culture cell receptor sites thereby blocking infectious particles. This can be a problem at low dilutions of tissue or ovarian fluid, particularly with the North American strain of VHSV.

**FBS** - fetal bovine serum taken from unborn calves in utero.

**Monoclonal Antibody (MAb)** - antibody produced by tissue culture cell lines derived from the spleen lymphocytes of immunized mice that have been fused (hybridoma) with mouse myeloma tumor cells. Hybridoma cells are cloned to select specific populations of cells, each producing a single antibody against one epitope or antigenic determinant site on one antigen molecule among those used to immunize the mice.

**Open System** - a system of incubating tissue culture cells that is open to the transfer of air, i.e., a plate. Requires a medium that is buffered against rising pH from air exchange. Common buffering systems are TRIS and HEPES.

**Overlay** - a medium used in the plaque assay that is placed over a virus-inoculated cell monolayer to prevent physical spreading of viral particles except by cell-to-cell release of infectious particles. The overlay contains a semisolid medium such as methylcellulose or gum tragacanth.

**Plaque** - a hole or focus of degenerate or dead tissue culture cells in the cell monolayer caused by viral replication. One discrete plaque is assumed to be caused by infection with one infectious particle or aggregate (called one plaque-forming unit = pfu).

**Polyclonal Antibody** - the entire population of antibodies produced in the sera of immunized animals that are directed against many epitopes on many of the antigenic molecules used for immunization. Most immunogens injected are whole cells or viruses that are composed of many different antigen molecules. Each antigen molecule may have more than one epitope. See "Monoclonal Antibody."

**Serum neutralization** - antibody molecules in the antiserum neutralize or block the antigenic receptor sites or otherwise degrade the protein coat (capsid) on the corresponding virus (antigen). This prevents virus attachment to and subsequent penetration of host tissue culture cells or virus replication once inside the cell. Neutralization of viruses by antibodies is specific and used to confirm viral identity. Neutralization may be reversible.

**Subculture** - transfer of inoculated tissue culture cells and supernatant from one plate to another that contains fresh cells. Used for suspected positive cultures to confirm presence of viral CPE as opposed to toxicity or contamination. Also used to replicate more viruses for storage, etc.

**TCID<sub>50</sub>** - denotes fifty percent tissue culture infective dose. This is the reciprocal of the highest dilution of virus that causes CPE in 50% of the wells inoculated with that dilution of infectious materials. This is determined by the Reed and Muench (1938) method.

**Tissue Culture-Grade Water** - High quality water (low in ions, minerals and contaminants) that must be used in preparation of all tissue culture media and reagents and in rinsing glassware to avoid toxicity to the cells.

**Titer** - the number of infectious units or plaque-forming units per unit of sample, i.e., per gram or ml.

**Toxicity** - changes in cell morphology or metabolism caused by toxic substances in the medium or inoculum. This can either cause cell death or interfere with cell metabolism, thereby reducing or preventing replication of the virus. These effects may have arisen through sample toxicity, bacterial or fungal contamination, improper glassware cleaning or improper media preparation. Usually toxicity can be distinguished from viral CPE by how rapidly it occurs (1 day), abnormal cell appearance without cell death, absence of the typical pattern of CPE for the test virus and, in the case of contamination, turbidity of the medium or visible contaminant colonies.

**NOTE:** Inoculation of very high-titer suspensions of certain viruses can cause an apparent toxic effect within 24 hours. If there is any doubt to whether disruption of the cell layer was caused by

toxicity or CPE, a subculture should be made. This is especially true for some inocula that can produce toxic effects that may take 5-7 days for development.

**Triturating** - The act of dispersing tissue culture cells for transfer by repeatedly drawing the cell suspension into a pipet and expelling it back into the flask. This should be done until the cells are in clumps of no more than three when examined with an inverted light microscope.

**Trypsin** - a proteolytic enzyme used to disperse cells and causes their release from the culture surface. Serum proteins neutralize it and its action is slowed by low temperature. Trypsin will cause release of the cells more readily than versene.

**Versene (EDTA)** - ethylene di-amine tetra-acetic acid is a chelating agent that binds divalent cations active in forming cell cement (hyaluronic acid) causing cells to round and release from the culture surface.



## Appendix B - Media Used in Tissue Culture & Virology

**NOTE:** All chemicals should be reagent or tissue culture grade. Use only glassware which is dedicated to tissue culture, or is new or has been acid washed.

### A. Antibiotic Incubation Mix

Anti-bacterial, anti-mycotic solution for disinfection of biological samples prior to inoculation of cell cultures in a plaque assay.

		<b>Solution Concentration</b>	<b>Final Concentration in Sample (Mixed 1:1)</b>
MEM-O	170 mls		
Penicillin/Streptomycin (10,000 units/ml penicillin, 10,000 ug/ml streptomycin)	40 mls	1600 IU/ml	800 IU/ml
Fungizone (250 ug/ml Amphotericin)	40 mls	40 ug/ml	20 ug/ml
Total Volume	250 mls		

Aseptically prepare the solution in 250 ml bottle.

Because the Fungizone is a saturated solution you must continually mix the cocktail mixture while dispensing. Dispense 1.0 ml volumes into 12-75 mm polypropylene tubes.

Snap caps closed. Freeze tubes for use as needed. Can store frozen for approximately 3 months. Avoid freeze-thaw cycles, thaw tubes immediately prior to use.

### B. Hanks Balanced Salt Solution (HBSS)

10x HBSS	100 mls
diH <sub>2</sub> O	895.3 mls
NaHCO <sub>3</sub> (7.5%)	4.7 mls

Mix and filter (0.2 µm) into small bottles

**C.** **MEM-0**

2.5x MEM	400 mls
L-glutamine (200 mM)	10 mls
NaHCO <sub>3</sub> (7.5%)	30 mls
Pen/Strep	10 mls
diH <sub>2</sub> O	450 mls
Tryptose Phosphate Broth	100 mls

Mix aseptically and filter (0.2 µm).

**D.** **2.5x MEM**

MEM (Eagle, modified)	9.4 g
diH <sub>2</sub> O	420 mls

Dissolve and autoclave 15 minutes at 121°C

**E.                                    MEM-10 / Hepes (Hepes Buffered Media)**

2.5x MEM	400 mls
Fetal Bovine Serum	100 mls
L-glutamine (200 mM)	10 mls
Hepes buffer (1 M)	15 mls
NaHCO <sub>3</sub> (7.5%)	5-10 mls
NaOH (1 M)	5 mls
diH <sub>2</sub> O	460 mls

Verify pH 7.3 - 7.6; filter 0.2 µm into bottles

**F.                                    MEM-5 Liquid Overlay (MEM-5/Hepes)**

2.5x MEM	400 mls
Fetal Bovine Serum	50 mls
L-glutamine (200 mM)	10 mls
Hepes buffer (1 M)	15 mls
NaHCO <sub>3</sub> (7.5%)	15 mls
NaOH (1 M)	5 mls
Fungizone	10 mls
Gentamicin	4 mls
diH <sub>2</sub> O	491 mls

Mix aseptically and dispense into bottles

### **G.                      Methyl cellulose (base for Overlay)**

Methyl cellulose (4000 centipoises)	8 g
diH <sub>2</sub> O	555 mls

1. Heat 225 mls diH<sub>2</sub>O in a 1 L bottle with a stir bar to a near boil.
2. Add methyl cellulose and mix by swirling.
3. Mix on a stir plate and add 330 mls cold diH<sub>2</sub>O washing down sides; stir until cool.
4. Stir at 4 °C overnight.
5. Autoclave for 15 minutes at 121 °C; will form an opaque solid.
6. Cool to room temperature and stir at 4 °C until soluble; store at 4 °C

### **H.                      Methyl cellulose Overlay**

2.5x MEM	400 mls
Fetal Bovine Serum	50 mls
L-glutamine (200 mM)	10 mls
Hepes buffer (1 M)	15 mls
NaHCO <sub>3</sub> (7.5%)	10 mls
NaOH (1 M)	5 mls
Fungizone (250ug/ml)	10 mls
Gentamicin (50mg/ml)	4 mls
Methyl cellulose	496 mls

Mix aseptically and dispense into bottles

## I.

### Trypsin-EDTA or TV

1. Trypsin-EDTA  
(10X), lyophilized  
0.5 g/L of Trypsin (1:250) and 8.5 g/L of NaCl

Rehydration: Rehydrate with 20 ml of sterile, distilled water.

2. Earle's Balanced Salt Solution (EBSS) (1X), liquid
3. Prepare a (1X) solution using EBSS

<u>Trypsin (1 part)</u>		<u>EBSS (9 parts)</u>
20 ml	+	180 ml
or 40 ml	+	360 ml
or 60 ml	+	540 ml
or 120 ml		1,080 ml

Take six 20-ml bottles (Trypsin-EDTA) and rehydrate with 20 ml/bottle of sterile tissue culture-grade water. Aseptically mix the 120 ml trypsin-EDTA with 1,080 ml of EBSS. Aliquot into 250-ml bottles. Label T/V Working Solution and freeze at -20°C.

An alternative procedure is: make versene-PBS using method described in K above. Rehydrate 1 vial of 2.5% Trypsin with 20 ml of tissue culture-grade water. After cooling versene-PBS, aseptically add 4 ml Trypsin into each of 5 bottles of 100 ml sterile versene-PBS. Final Trypsin concentration is 0.1 %